



LiveLife Watch

Your
mobile alarm is
ready to use



- ✓ SIM card installed
- ✓ Alarm programmed
- ✓ Alarm is ready to use

USER MANUAL

Stay Safe, LiveLife

www.livelifealarms.co.nz

0800 555 179



Accessories included

Magnetic charger base



USB charging cable




USB/240V power adapter



Bluetooth Beacon



Getting to know your LiveLife Watch




Please note the device can support multiple beacons if required.

*Watch comes with SIM card installed.



Charging your LiveLife Watch

- Connect the USB cable from the charging base to the 240 Volt AC adaptor.
- Place the watch into the magnetic charging base correctly. You will hear an announcement: *Your LiveLife Alarm is charging.*
- When the battery drops to 20% charge an announcement from the device will occur advising low battery.
- Upon an announcement, please fully charge your battery.



LiveLife
Watch

SOS**Activating an SOS alarm**

- Press the side button for 1–3 seconds until you feel your watch vibrate.




- When activated, you will hear an announcement confirming the SOS activation. You will have 10 seconds to cancel the alert. To cancel the activation, simply tap the cancel icon on the screen or click the watch side button. If you do not cancel the SOS activation, the alarm will immediately notify your elected emergency contacts.

Family and Friends Monitoring Option:

The watch will send a text message to the programmed contacts (up to 6) with a link to your location on Google Maps, then it will start the call sequence.

There is a 10 second delay between each contact, you can stop the call sequence during this time by tapping the cancel icon on the screen.

**Testing the SOS alarm**

To ensure your watch operates correctly, we recommend testing it as soon as you receive it. To perform the test, simply activate the SOS alarm by pressing the SOS button. Before doing so, notify your emergency contacts that you are conducting a test and provide them with your alarm's phone number. It is best practice to ensure the receipt of both text messages and call to emergency contacts when testing the alarm.

**Fall detection**

Your LiveLife Watch can automatically detect a fall. When this occurs you will hear an announcement confirming the fall alert activation. You will have 10 seconds to cancel the alert.



If you do not cancel the fall alert activation, the alarm will immediately notify your elected emergency contacts, according to your monitoring model.

Testing fall detection

When you first receive your LiveLife Watch, you can gently drop it or toss it onto a soft surface, such as carpet, to observe the fall alert sequence in action. However, we advise against testing the fall detection feature repeatedly. Each fall is unique due to various factors, and excessive testing may result in injury or damage to the device. Never intentionally drop or throw the watch onto a hard surface to test this feature.

What can stop fall detection from working?

The fall detection feature uses a 3D accelerometer to measure the speed of impact, which means falls from certain heights may not trigger an alert. For instance, a fall from less than one metre, such as from a seated position, may not meet the criteria for a heavy fall. Additionally, false alarms can occur. If needed, you can adjust the sensitivity settings or turn off the fall detection feature by calling our support team on **0800 555 179**.



Never intentionally drop or throw the watch onto a hard surface to test this feature.

**Receiving calls on the LiveLife Watch**

Anyone who knows the device's phone number can call it. To answer an incoming call, tap on the green icon on the watch screen. You can find the device's phone number on the packing slip included with your LiveLife Watch.


**Finding the wearer's location**

The watch uses Bluetooth, WIFI and GPS to determine its location. You can send the text command **LOC** to the device's phone number and it will reply with its location on Google Maps.



Other watch functions

When the watch is turned on, slide left or right on the watch face to view available functions including Steps, Heart Rate, Contacts and Settings. Press on the icon to view the options for each group of related settings. The device will provide an audible confirmation when an option is selected.



Time & Date

Time can be displayed in analog or digital format.

Step Counter

Tracks Daily Steps, Distance, and Kilojoules Burned*

Heart Rate

Heart Rate Monitor: Measures heart rate in real-time*.

Contacts

A call can be made to an emergency contact within the limits of our Emergency Use Policy**.

Settings

Options to change the watch faces, turn off and restart the watch.

*Measurements are for wellness purposes only and not intended for medical use.

** Available on page 6 of this document.



Maintaining your device

Using your watch in wet areas

Your watch can be used in the shower and other wet areas. Whilst it can be immersed briefly, it is not intended for long-term immersion and especially not in deep water or water which is significantly hotter or colder than the ambient temperature. The product should be dried thoroughly after use in wet areas. Do not charge the product until fully dry.

Cleaning Recommendations

Use only mild soap and soft cloths to clean or dry your watch, and rinse with water afterwards. Do not try to insert anything into the speaker or microphone apertures. Do not use any solvents, especially isopropyl alcohol to clean the glass. Do not use pressure cleaning (jets etc.) or jewellery cleaning ultrasonic baths in particular, regardless of whether cleaning chemicals are included as the mechanical action can cause damage.



Monthly alarm health check

Performing regular tests of the LiveLife Watch

Like other electronic devices, it is important to conduct regular checks of the watch to make sure it is working correctly. We suggest a monthly health check which includes:

- Ask a relative or friend to call the alarm's phone number and you can answer in hands-free, two way voice.
- Check your watch for obvious signs of damage. Make sure the four brass discs on the back of the watch are bright and shiny.
- Activate the SOS then wait until you hear the announcement. Click the SOS icon on the watch face to cancel the activation or press the watch side button.



Mobile network & SIM card

Ready to go out of the box.

As part of the setup and pre-programming we have activated and installed your nano-sized SIM card. Your device is ready to go!

Changing or replacing a SIM card

The SIM card should only be removed by LiveLife Alarms or under instructions from us.



After sales support

Free lifetime support

Your watch comes with free lifetime support. If you have any concerns please call **0800 555 179** to speak with a customer support consultant. Our customer support opening hours can be found on our website at:
www.livelifealarms.co.nz/contact-us/



Specifications

Dimensions: 47.5mm x 16.9mm x 1.87mm.

Weight: 65g.

Battery Type: Rechargeable, 3.7V, 580mAh.

Battery Duration: Up to 36 hours under typical use.

Charging Voltage: 5V DC.

Showerproof Rating: IP67.

Locating Technologies: Bluetooth, WiFi, GPS.

Screen Type: OLED 0.96 inches.

Eyeglass: Toughened organic glass.

SIM Card: Spark

NBN, Internet, and WIFI:

The LiveLife Watch does not require NBN or a landline to function and does not need an internet connection. It may use WIFI for location purposes but is not reliant on WIFI.



How to make changes to your emergency contacts

Available on family and friends monitoring plan only.

You can make changes to the emergency contacts and other settings by sending simple commands to the alarm's phone number via text from a mobile phone. To change an emergency contact number send these commands. Do not include spaces in any text command. Is not case sensitive. The watch will respond to each text command with a reply text to confirm the change.

Person 1		A1,1,1,0400123456,Name	Send
Person 2		A2,1,1,0400123456,Name	Send
Person 3		A3,1,1,0400123456,Name	Send
Person 4		A4,1,1,0400123456,Name	Send
Person 5		A5,1,1,0400123456,Name	Send
Person 6		A6,1,1,0400123456,Name	Send



How to change the watch face




1. From your current watch face, swipe left until you reach the Settings screen.
2. Tap the screen to enter the Settings menu.
3. Select "Watch Faces" to open the watch face options.
4. Swipe left or right to browse through available designs.
5. To see more variations of a watch face, tap the gear (cog) icon at the bottom of the selected design and swipe up or down to see the additional options.
6. Tap your desired watch face to apply it.



Getting to know Bluetooth Beacon


Your Bluetooth beacon is already turned on and ready to use.

Top



Hold down the ON/Off Button for 3 Seconds to turn on or off.

Back



Area to stick double sided tape, to make it harder to lose your beacon.

Side



Your Bluetooth Beacon

The LiveLifeWatch comes with a Bluetooth beacon to place inside your home. The beacon has been programmed by us with the exact co-ordinates of your home address you supplied with your order. When you activate the SOS on your LiveLifeWatch while at home it will look for the beacon and those programmed co-ordinates so it can send that exact location in its help text messages via Google Maps.

If you are outdoors, wearing your LiveLife Watch and out of range of Bluetooth it will use GPS, or Wifi to get your location.

Leave your beacon at home

Your beacon is meant to be left at home. If you move house please call us on 0800 555 179 so we can arrange to re-program your beacon. If you holiday in another location purchasing a second beacon could be a sensible option.

Testing the accuracy of your beacon

You can test your beacon by sending a SMS to your LiveLife Watch from a mobile phone. Send the SMS command: **LOC (not case sensitive)** to the device's phone number.

If the Watch and beacon are on and the beacon is within range of the device the LOC response will state:
Location by Beacon

Plus also your location and time information.

Features

1. Compact and lightweight design for easy installation.
2. Designed for use with LiveLife personal mobile alarms.
3. Includes a C2477 button battery, with an estimated lifespan of approximately two years.
4. Beacons offer both home and away information, enhancing the functionality of the mobile alarm systems.



Important Information

Emergency usage policy & excessive use

The LiveLife Watch is designed to be used in an emergency and not as a communication device such as a mobile phone. Your watch comes under our 'Fair Use Policy' which means it is to be used for help in emergency situations. We have allocated a fair and adequate amount of voice and data to your device based on average usage patterns. Voice and text usage for every device is monitored by LiveLife Alarms. Customers using the devices contacts to make outgoing calls heavily will be contacted to discuss usage options. We reserve the right to apply further charges for excess use or discuss the option to move the SIM card into your name.

Spark Mobile connection

Your watch is relying on a connection to the Spark 4G mobile network to perform all of its emergency text and call functions. It also uses Bluetooth, GPS and WIFI to perform its location functions.

Spark Mobile SIM card activation

We have installed a Spark SIM card and activated it. Please refer to your packing slip for the details of the programmed emergency contact numbers.

12 months warranty

The LiveLifeWatch comes with a standard 12 months warranty from date of purchase. If the watch has a manufacturing defect causing it to fail within the first 12 months of purchase, we will replace it free of charge or refund the purchase price.

How to claim against your warranty

Please contact us by emailing us at: info@livelifealarms.com.au or call us on 0800 555 179 to arrange to have the watch returned for a warranty claim assessment.

Circumstances that can void your warranty

Please be aware your warranty can become void if any of the following have occurred:

- The watch has been opened or tampered with.
- The watch has been modified or not used in accordance with our instructions.
- The watch has been dropped onto a hard surface or suffered damage from physical impact.
- The watch has been damaged by fire or has been exposed to excessive heat.
- The watch has been used with an incorrect power adapter.
- The watch has been subject to an incorrect voltage due to lightning or a power surge.
- The watch has been subject to a liquid spillage causing a failure.
- The watch has been submerged in water.

What if something goes wrong after 12 months?

If the watch breaks after the first 12 months please contact LIFELIFE ALARMS to send it back to us for repair or replacement. We will only charge you for the cost of repair and shipping. You won't have to pay full price for a watch.

Guarantee under Australian consumer law

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

LIFELIFE ALARMS is not liable for any indirect or consequential loss or claims.

Refunds & Returns Policy

30 day refund and return policy

If for any reason you are not completely satisfied with your purchase we will give you a refund from the date you receive the goods till the end of 30 days from that date.

Please email us at: support@livelifealarms.com.au within that time or call us on 0800 555 179 if you are not satisfied with your purchase so that we can resolve any problems.

This refund policy does not apply to goods which show wear and tear or are damaged after delivery or if any attempt has been made to alter the product or if they have been dropped or broken. All products and parts must be returned in their original condition. All postage and insurance costs are to be paid by the user. We recommend that you return the product via Registered post and that you pre pay all postage. The user bears any risk of lost, theft or damaged goods during transit and therefore we advise you take out shipment registration of insurance with your postal carrier. LIFELIFE ALARMS will not be responsible for parcels lost or damaged in transit if you choose not to insure.

Refunds (Coverage)

It is ultimately the user's responsibility to test their alarm during the 30 day refund period and identify any mobile network coverage or signal issues.

We use our best endeavors to ascertain if there is adequate coverage at the location where the user resides. Our standard practice before setting up your watch is to do a coverage check on the Spark Mobile Network' at the supplied, primary address where the watch will be used if and if necessary also ask the purchaser if they personally are aware of any weak signals inside the residence on devices connected to the Spark Mobile Network.

Like any mobile device, the LiveLife Watch relies on mobile network coverage, which may vary by location and building structure. While we conduct a coverage check prior to shipping, we cannot guarantee uninterrupted service in all environments, especially indoors or in remote areas. If after testing the user believes there is a weak signal then they should contact us via phone or email before the '30 day period' expires in order for us to investigate the coverage issue.

If after investigating we find the coverage is inadequate and we cannot provide an alternate SIM card and network provider, upon the device being returned we will refund the purchase price.

Medical & legal disclaimer

General information about medical matters

The LiveLifeWatch 'User Manual' may contain general information about medical and health matters. The information is not advice, and should not be treated as such.

Limitation of warranties

The medical and health information within this 'manual' is provided "as is" without any representations or warranties, express or implied. LIFELIFE ALARMS and its employees and/or agents make no representations or warranties in relation to the medical and health information within this manual.

• Without prejudice to the generality of the foregoing paragraph, LIFELIFE ALARMS does not warrant that:

- The medical and health information on its website or in this 'manual' will be constantly available, or available at all; or
- Any medical and health information on its website and within this 'manual' is complete, up to date or error free to the best of our knowledge.

Liability

You must not rely on the information within this 'manual' as an alternative to medical advice from your doctor or other health services provider.

If you have any specific questions about any emergency plan, medical or health matters you should consult your doctor or other professional health services provider or personal carer. You should never delay seeking medical advice, disregard medical advice, or commence or discontinue any medical treatment because of information within this 'manual'.

Liability limitations

Nothing in this legal disclaimer will limit any of our liabilities in any way that is not permitted under applicable law, or exclude any of our liabilities that may not be excluded under an applicable law.

Product Disposal

Your LiveLife personal alarm product contains an embedded lithium-ion battery. These should not be disposed of as general waste or recycling, as the battery could present a fire hazard if pierced or damaged, and the internal PCBA contains recyclable elements. Depending on state, various options are available to you. We recommend consulting the guidelines offered by your state, such as the NSW EPA, for community recycling stations that accept these products.

Transporting or shipping your product by air

If you need to ship your product by air, we recommend allowing the battery state of charge to fall below 40% before you securely pack the product for shipping.



LiveLife
Watch

www.livelifealarms.co.nz

0800 555 179